

AY-A2000

DENTAL CHAIR UNIT



Effective Treatments Are At Your Fingertips

To create a comfortable clinic environment

- Ergonomic Design
- More Stable And Sturdier
- Comprehensive Treatment



ANYA A2 CONTROL SYSTEM

Multifunction

Emergency position

one-step spitton position

nine memory positions

X-RAY film viwer



Bringing A comfortable Experience During The Operation



New motor control program

optimized control algorithm

flexible motor speed changing

To bring a more silky treatment
experience to the patients



ERGONOMIC DESIGN

CONSIDERATE DESIGN

Butterfly shape backrest , fit the curve of the back;
Lumbar uplift, multi-angle adjustment of the head pillow, supporting the waist and neck

High elastic inner core
Powerful soft elastic

Skin friendly fabric
delicate touch

Upgrade process
anti-bacterial and anti-mildew



Multifunction Assistant Tray

Which is convenient for the assistant to assist the dentist for comfortable treatment



New Design Fashion Cabinet

Antioxidant design



More Stable And Comfortable

Thickened base, butterfly shape metal backrest equipped with air spring, which increase the stability of the dental chair



Nebula-FS

Luxury inductive LED cold lamp
can adjust brightness

5 Color
Temperature Modes



White Light Mode



Yellow Light Mode



Mixed Light Mode

"Nebula" Series Reflective Oral Lamp

Inductive switch

- Brightness 7000-42000 Lux, stepless adjustment
- Different color temperature mode: white light, yellow light, mixed light,
- Precise digital display of brightness and color temperature
- Color rendering index > 90



Selected quality materials

High configuration
to create standard quality



Medical Grade Durable
Corrosion-resistant Water Pipe



Taiwan Timotion Motor



ODE Italy Imported Solenoid Valve



Multifunction Instrument Tray

Equipped with rotatable tray

Rotatable Spittoon

Optional glass/ceramic spittoon
Easy dismounting, easy cleaning





Low Pedal, Easy To Step On Work In High Efficiency

The Transmission Structure Is Optimized,
And The Trampling Is More Comfortable, Less Noise
Lower Failure Rate, Longer Service Life