

安雅医疗  
Anya Medical

# AY-A3000

## DENTAL CHAIR UNIT





# Effective Treatments Are At Your Fingertipst

- Ergonomic Design
- More Stable And Sturdier
- Comprehensive Treatment





# ANYA A2 CONTROL SYSTEM

## Multifunction

Emergency position

one-step spitton position

three memory positions

X-RAY film viwer





# Bringing A comfortable Experience During The Operation



New motor control program  
optimized control algorithm  
flexible motor speed changing  
To bring a more silky treatment  
experience to the patients





# ERGONOMIC DESIGN

## CONSIDERATE DESIGN

The design of the seat is based on the shape of the human body to provide the patient with a sense of comfort and support

Comfortable design of The chair cushion

Environmental protection material





# Multifunction Assistant Tray

Which is convenient for the assistant to assist the dentist for comfortable treatment





**Flexible adjustment  
of chair position**

Easily adjust for  
better position.



# Large Capacity Tissue Box

Intimate design, which is more comfortable for patients





# Nebula-FS

Luxury inductive LED cold lamp  
can adjust brightness

5 Color  
Temperature Modes



White Light Mode



Yellow Light Mode



Mixed Light Mode



# "Nebula" Series Reflective Oral Lamp

Intelligent touchless control

- Brightness 7000-42000 Lux, stepless adjustment
- Different color temperature mode: white light, yellow light, mixed light,
- Precise digital display of brightness and color temperature
- Color rendering index > 90





# Selected Quality Materials

High configuration  
to create standard quality



Medical Grade Durable  
Corrosion-resistant Water Pipe



Taiwan Timotion Motor



ODE Italy Imported Solenoid Valve





# Adjustable Headrest

Supports the back of the neck



# Sturdy medal backrest

Can bring you a relaxing and comfortable treatment feeling







# Low Pedal, Easy To Step On Work In High Efficiency

The Transmission Structure Is Optimized,  
And The Trampling Is More Comfortable, Less Noise  
Lower Failure Rate, Longer Service Life