

Steam Sterilizer Autoclave

Class B | 23Litre



Smart Operation



Smart Operation





Large Capacity. Open stainless steel water-tank
Easy clean No deformation



Based on the new class B of international prEN13060
have 3 pulsating vacuum, Vacuum up to -0.9Bar (-0.09Mpa)
Suitable to sterilizing wrap, unwrap, solid, porous, hollow device.



Large Capacity Liner.304 Stainless Steel Tray



with built-in printer .USB optional
The disinfection process data and date automatically print.



Use imported&advanced16-bit microprocessor,easy to operate.
Most suitable for the Dept of Stomatology,Ophthalmology,
Surgery and Lab.



Double lock door system for safety during operating.